
12/96-Hu/Brs-

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The I-measuring amplifier D is a DC amplifier; when used in conjunction with a voltmeter as the display unit, it serves to measure currents, voltages and charges.

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#### 1 Safety notes

- To use the device with a mains voltage of 115 V AC, convert it as described in section 4.2.
- Use only with non-hazardous contact voltages.
- Always start with the measuring-range selector ⑧ switched to the most insensitive range, particularly when the magnitude of the measured value is unknown.
- Hum signals above 10 V<sub>p</sub> at output ⑧ falsify the measuring result. You may need to check this (e.g. with an oscilloscope)!
- Input currents greater than approx. 5 mA (e.m.f.  $\geq 500$  V) overdrive the amplifier. Use a resistor if necessary! (see section 3.4).
- Always observe the instructions in sections 3.2 to 3.3 to avoid falsified measurements.

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## 2 Scope of supply, description, technical data

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- ① Mains switch with on-indicator lamp
- ② BNC input for voltage measurements
- ③ BNC input for current and charge measurements
- ④ Ground socket of input, connected to protective-earth conductor
- ⑤ Zero-key. The input of the measuring amplifier is grounded
- ⑥ Pair of output sockets for connecting a voltmeter measuring range 0.3 V DC to 10 V DC; see table 2)
- ⑦ Zero-point adjuster

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- ⑧ Measuring range selector
- ⑨ Screened cable BNC/4-mm (575 24) with separate connector (4-mm plug) for screen  
Impedance: 50 Ω  
Capacitance: 120 pF  
Length: 1.2 m

## 2.2 Technical data

Current ranges	Input resistance	Voltage ranges	Input resistance	Charge ranges	Input resistance
10 <sup>-6</sup> A	10 <sup>5</sup> Ω	10 <sup>-2</sup> V	10 <sup>9</sup> Ω	10 <sup>-6</sup> As	10 <sup>5</sup> Ω
10 <sup>-7</sup> A	10 <sup>5</sup> Ω	10 <sup>2</sup> V	10 <sup>9</sup> Ω	10 <sup>-7</sup> As	10 <sup>5</sup> Ω
10 <sup>-8</sup> A	10 <sup>5</sup> Ω	10 <sup>1</sup> V	10 <sup>9</sup> Ω	10 <sup>-8</sup> As	10 <sup>5</sup> Ω
10 <sup>-9</sup> A	10 <sup>5</sup> Ω	10 <sup>0</sup> V	10 <sup>9</sup> Ω	10 <sup>-9</sup> As	10 <sup>5</sup> Ω
10 <sup>-10</sup> A	10 <sup>5</sup> Ω	10 <sup>-1</sup> V	10 <sup>9</sup> Ω		
10 <sup>-11</sup> A	10 <sup>5</sup> Ω				

Table 1

All data refer to an output voltage of 1 V. However, the device can be overdriven up to an output voltage of 10 V<sub>p</sub>.

**Input:** one BNC socket each for current/charge measurement and voltage measurement. Ground electrically connected to PE conductor via PC board.

**Output:** two 4-mm sockets, source resistance 100 Ω  
**Measuring accuracy:** typical 3 % for current/charge measurement  
 5 % for voltage measurement  
**Zero-point stability:** 0.1 %  
**Gain stability:** error < 1 %

**Mains connection voltage:** 230 V 50/60 Hz, convertible to 115 V AC (see section 4.2)  
**Power consumption:** 6 VA  
**Fuses:** T 0.125 B for 230 V AC  
 T 0.2 B For 115 V AC  
**Dimensions:** 20 cm x 21 cm x 23 cm  
**Weight:** 2 kg

## 3 Operation

### 3.1 Putting into operation, selecting the display instrument

Connect a voltmeter at output ⑥. All DC measuring instruments with an internal resistance greater than 3 kΩ and at least one measuring range between 0.3 V and 3 V (maximum 10 V) can be used, e.g.:

Device	Cat. No.	Measuring ranges			
		0.3 V	1 V	3 V	10 V
Moving-coil instrument P	531 80		X	X	
Demo-multimeter	531 911 or 915	X	X	X	X
AV-meter	531 94		X	X	X
Oscilloscope, e.g.:	575 21	X	X	X	X
Xy-yt-recorder, e.g.:	575 663		X		X
Yt-recorder, e.g.:	575 702		X		X
for computer-assisted experimenting: e.g. CASSY-E	524 007	X	X	X	X

Table 2

Inadequate screening (due to the experiment conditions) can result in a 50 Hz signal being superpositioned on the quantity to be measured; this signal is also amplified, even though its display is not desired.

The above measuring instruments average a symmetrical hum to approx. zero (via a parallel-connected capacitor if necessary). The hum signal at the output may not exceed 10 V<sub>p</sub>, as otherwise it is asymmetrically limited and thus causes

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an erroneous reading. To check this, you may need to connect an oscilloscope at output ⑥.

### 3.2 Current measurement

Connect the screened cable at BNC input ③ and supply the measurement signal.

Set the current measuring range on ⑧ so that (depending on the measuring instrument) the output signal connected is about 0.3 V to 3 V (maximum 10 V).

The measured value is determined by multiplying the measuring range by the numerical value of the output voltage (measured in volts).

Example:	Measuring range:	$10^{-8}$ A
	Output voltage:	1.3 V
	Measured value:	$1.3 \cdot 10^{-8}$ A

#### Important:

- Screen the input as well as possible to avoid input interference, e.g. from the mains frequency. When conducting stationary current measurement, a capacitor, e.g. 0.5 µF (538 04), can be connected parallel to the screened-cable input for smoothing (see also 3.1).
- Do not disturb the screened cable during the experiment to prevent charge shifts.
- Isolate the measuring point from other current sources (e.g. also from the insulation of voltage-bearing leads and the benchtop) to minimize leakage currents. Compensate any remaining input offset currents using the zero-point adjuster on ⑦.

### 3.3 Voltage measurement

*Attention:* use only non-hazardous contact voltages (see section 1)

Connect the measuring lead at BNC input ② and supply the measurement signal.

Set the voltage-measuring range on ⑧ so that (depending on the measuring instrument) the output signal connected is about 0.3 V to 3 V (maximum 10 V).

The measured value is determined by multiplying the measuring range by the numerical value of the output voltage (measured in volts).

Example:	Measuring range:	$10^1$ V
	Output voltage:	0.7 V
	Measured value:	7 V

### 3.4 Charge measurement

#### Important:

The input current may not be greater than approx. 5 mA, as otherwise the amplifier is overdriven and the measuring result is falsified. Therefore, **at e.m.f. values greater than 500 V** an additional resistor of e.g. 1 MΩ (for voltages up to 5 kV) must be connected in front of the amplifier input.

Connect the measuring lead at BNC input ③ and supply the measurement signal; for e.m.f. values over 500 V, always use an additional resistor of e.g. 1MΩ.

Set the measuring range so that (depending on the measuring instrument) the output signal connected is about 0.3 V to 3 V (maximum 10 V).

The measured value of the charge is determined by multiplying the selected measuring range by the numerical value of the output voltage (measured in volts).

Example:	Measuring range:	$10^{-8}$ As
	Output voltage:	3 V
	Measured value:	$3 \cdot 10^{-8}$ As

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
Fig. 2.1


Fig. 2.2

## 4 Replacing the fuse, converting the mains voltage

### 5 Replacing the primary fuse

Pry out insert (a) with holder for primary fuse (b) and spare fuse (c) (Fig. 2.1).

Replace the defective fuse (b) using a new fuse (c) which has been checked for the correct rating (Fig. 2.2).

Insert the new spare fuse (c) (see section 2 for the correct rating) and put the insert (a) back in the device.

### 5.1 Converting the device for operation with 115 V mains voltage

*Attention: always unplug the device first!*

Loosen housing screw (e) on the bottom of the device using a Philips screwdriver (size 2) (Fig. 3.1).

Turn the device right side up and remove the top housing section (f) (Fig. 3.2).

Remove the PC board (h) in the 230 V position from the blue slot (g) at the transformer (Fig. 3.3).

Turn the board (h) around so that the stamped designation "110 V" (for 115 V mains voltage) is at the bottom left above the "1" on the contact side of the plug-in slot (g) (Fig. 3.4).

Insert board (h) and screw the housing back together.

Insert the correct primary fuse for the new mains voltage (see technical data).

See section 4.1 (Fig. 2.1 and 2.2) for instructions on how to change the fuse.


Fig. 3. 1 - 3.4