3/92 -Sf-


Fig. 1 Motor zum Röntgengerät 42 kV

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Fig. 1 Motor for the X-ray apparatus 42 kV

The motor is used for driving the counter tube goniometer in the X-ray apparatus, 42 kV (554 90). By outputting an angularly proportional voltage in conjunction with a counter tube apparatus, the automatic recording of the X-ray spectra is permitted. The spectra are either plotted with a recorder or displayed on the monitor (and, if required, output on a printer) of a computer measurement station equipped with an interface. The control of the motor (speed and rotation direction) is carried out internally or by means of an external voltage, which is supplied, for example, by an interface which is program-controlled.

7 materials of various atomic number ($23 \leq Z \leq 49$) are contained in the set of absorption films (554 925) for the recording of absorption spectra (required, e.g. for the investigation of Moseley's law and the λ^3 -law for atomic attenuation coefficient).

1 Safety instructions

- As no mechanical alterations are carried out on the housing, the design of the X-ray apparatus 42 kV (554 90) is not infringed upon when the motor is assembled according to the prescribed stipulations. Consequently, the license included with the X-ray apparatus remains valid.
- Maximum permissible voltage supply at the input ④: 14 V AC
- Maximum permissible control voltage at the input ③: ± 8 V DC
- Directly before initial assembly, remove protective foil from the adhesive surface of the velcro band ⑦.
- Handle the absorption films carefully.

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2 Description; technical data; scope of supply

2.1 Motor for X-ray apparatus kV (554 905); see Fig. 1

Precision motor with coupled 10-turn precision potentiometer (10 k Ω) connected via slide-clutch coupling, transmission ratio (1:834) for the drive adaptor ① and control electronics built-in to housing

Motor data:

Control voltage:	max. 8 V
Torque:	max. 60 N cm
Current consumption at nominal voltage and nominal load:	330 mA

- ① Selection switch for internal and external motor control
Internal control: positions L and LL for slow and fast counter-clockwise rotation
positions R and RR for slow and fast clockwise rotation
Speed change for slow operating mode with controller ②
external control via voltage at input ③: Control ext.
- ② Speed controller for internally controlled slow, counter-clockwise and clockwise rotation (positions L and R of switch ①)
- ③ 4-mm input sockets for external control voltages (max. ± 8 V) (for position "ext." of switch ①);
Lower socket connected to ground
- ④ Jack bush for cable jack ⑩ for the connection of the supply voltage
Required voltage: 12 V AC (max. 14 V AC)
- ⑤ Output sockets for angle proportional voltage (max. ± 10 V)
Lower socket connected to ground
for the angle range 0...30°/60° of the counter tube goniometer in the utilizable voltage range of the X-ray apparatus:
-3 V... +2.4 V
Smallest utilizable angle step: 0.1°
- ⑥ Drive adapter for attachment to the detector control knob on the X-ray apparatus: with washer (6.1) for transmitting the torque
- ⑦ Velcro band for mounting of the motor to the housing of the X-ray apparatus
- ⑧ Bore hole for the supporting rod ⑨
- ⑨ Supporting rod, to be screwed into the hole ⑧; with lock nut (9.1)
- ⑩ Two-core cable for connection of the supply voltage to the jack bush ④

2.2 Set of absorption films (554 925)

7 films, suitable for mounting in the film holder of the collimator in the X-ray apparatus 42 kV

Z	Element		Length	Thickness
13	Al	Aluminium	40 mm	0,5 mm
26	Fe	Iron	40 mm	0,5 mm
29	Cu	Copper	40 mm	0,07 mm
40	Zr	Zirconium	40 mm	0,05 mm
42	Mo	Molybdenum	33 mm	0,1 mm
47	Ag	Silver	33 mm	0,05 mm
49	In	Indium	33 mm	0,3 mm
			Width: 15 mm	

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3 Operation

Note:

Calibrate the motor before assembly.

Remove protective foil from the adhesive surface of the velcro band ⑦ before initial assembly!

If necessary, grease the washer (6.1) of the drive adapter ⑥ with a thin layer of Vaseline, or something similar, to ensure that it slides smoothly.

Pay careful attention to the information from the instruction sheet provided for the X-ray apparatus 42 kV (554 90).

3.1 Experiment arrangements

The spectra can be plotted either with an XY-recorder or computer, displayed on a monitor and output on a printer (Fig. 2.1/2).

In both cases the following equipment is required:

X-ray apparatus 42 kV (554 90) with X-ray tube (554 94) and monocrystal (554 77 or 554 78)

End window counter for beta, gamma, and X-rays (559 05)

Cable for counter tubes, 1 m (559 07)

AC power supply, 12 V (e.g. 562 73)

For plotting the spectra use either
the ratemeter with counter tube connection and analog output (e.g. 575 52)
XY recorder (e.g. 575 662)

or
MS-DOS computer (640 kB) with MS-DOS operating system, version 3.0 or higher,

Cassy components:

- MS-DOS connector S or L (524 001 or 002), installed in computer,
- CASSY-E or P (524 000 or 728 307)
- GM-box (524 033); alternatively: counter device with voltage supply for counter tube and pulse output (e.g. counter S 575 47, with plug-in power unit, e.g. 530 88)
- Software "Bragg" (524 842)

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Fig. 2.1 Computer-assisted recording of spectra

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
Fig. 2.2 Plotting spectra using an XY-recorder

Fig. 3 Attaching the motor to the X-ray apparatus
 Note: attaching half (c) of the velcro band by means of adhesive with correct positioning during initial assembly of the motor (see section 3.2)

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3.2 Calibration; assembly; recording measured values

Note: Information on calibration and recording measured values for computer-assisted recording, as specified in Fig. 2.1, is provided in the program "Bragg" as well as in the accompanying brochure.

Screw the supporting rod ⑨ into the threaded hole ⑥ as far as possible (for computer operation only a f t e r the program "Bragg" has been loaded) and connect the power supply (12 V AC) to socket ④.

In order to match the motor to the angle range (0°...30°/60°) of the counter tube goniometer, set the control voltage b e f o r e assembly so that the output voltage at socket pair ⑤ lies slightly above -3 V (approx. -2.9 V).

As specified in Fig. 3, attach the motor to the detector knob (a) of the coupled and carefully adjusted goniometer using the drive adapter ③ so that the velcro band runs directly along edge (b), i.e. the two halves (c) and (d) of the velcro band fit exactly one on top of the other with no overlap (after initial attachment, the velcro band half (c) remains attached to the X-ray apparatus due to adhesive).

Tightly press the motor housing in this position against the housing of the X-ray apparatus;
 Screw in the rod ⑨ just far enough so that the motor housing is supported without tilt and fasten with the lock nut (9.1).

To plot the spectrum using the XY-recorder (Fig. 2.2), simultaneously start the recorder and the motor (selector switch ① in position L); for computer-assisted recording follow the instructions given in the program "Bragg".